

**SINGLE USE**



# DIAGNOSTIC Catalogue

# KEEPING IT SIMPLE

- > PRECISE
- > QUALITY
- > ERGONOMIC
- > STERILE
- > SINGLE USE



# CONTENTS

- 4 ABOUT US
- 5 OUR PRODUCTS
- 6 WASTE MANAGEMENT
- 8 TWO PIECE  
DIAGNOSTIC KITS
- 12 THREE PIECE  
DIAGNOSTIC KITS
- 20 COMPLETE PERIODONTAL  
DIAGNOSTIC KIT
- 22 PACKAGING
- 23 MATERIALS
- 24 HOW TO BUY
- 25 SHARPS CONTAINERS

# ABOUT US

Our goal is to simplify dentistry as you know it. Minimising instrument preparation into one easy step, our single-use instruments will streamline the way you work. Quick, convenient and stress-free, they are improving dental clinics on a daily basis.

Pre-packed and sterile, they arrive ready for use. Removing the need for sterilisation allows you to focus on what's important: delivering high-quality patient care.

We know running an efficient clinic can be a juggling act. Our instruments will save your team time and take the stress out of a busy schedule.

Join us and take a convenient approach to dentistry.





# OUR PRODUCTS

By listening to customers and continuously working to improve our offering, our range has evolved to feature a wide variety of instruments and kits ideal for every procedure.

## Individual Instruments

Using individually wrapped instruments is a great way to transition your clinic towards implementing single-use solutions. Only open those you need for procedures and take minimising waste to a whole new level.

## Diagnostic Kits

Have you been entertaining the idea of single-use instruments? Equipping your clinic with our readily available diagnostic kits allows your nurse to spend less time setting up and more time chair side, focusing on the patient. Who would complain about that?

## Restorative Kits

Our restorative instrument kits will leave you staring at the autoclave wondering why you spend so much money maintaining it every year. Using our thoughtfully designed combinations alleviates the pressure and expense of reprocessing. Remember if you can't find a restorative combination that works for you, we'll make you one. That simple!



# WASTE MANAGEMENT

Our disposable instruments solve many problems for the dental industry. Reduced exposure to infection, minimised reprocessing and complete convenience... they really are the obvious choice.

Our mission is to ensure that our instruments become a single-use, zero waste solution. Through determination and collaboration with stakeholders who have a similar goal we have made significant progress in recent years.

The current waste management program is focused on a method that utilises a neutral waste recycling solution. The process minimises carbon emissions and instead repurposes the instruments into an energy source in the form of electricity, heat or fuel.

As we implement this proposal and develop a better understanding of our waste potential, we will be able to utilise this information to develop relationships with companies that want to recover and repurpose the raw material resources. This is our future vision and goal.




We are now offering a SMART Waste Recovery Program which enables clinics to dispose of their instrument waste into an allocated sharps container, have it collected, then taken to be processed within a dedicated waste facility.

To take advantage of our program all you need to do is register your interest via our website [here](#) or simply call our team on 07 5445 6486 and we can chat further to see if this program is suitable for your requirements. You can also find some more information [here](#).

The program includes the supply to clinics, collection when full and disposal of the waste.

See how the process works below.



\* Browse our sharps container range towards the back of this catalogue.

# TWO PIECE DIAGNOSTIC KITS

Our two-piece diagnostic kits are ideal for an initial check of the mouth. The anti-fog mirror delivers clear vision and our combined probes give you two instruments in one. Browse our selection to see which mirror and probe combination suits you best.





## TWO PIECE DIAGNOSTIC KITS



**D2004**

### Diagnostic Kit

Great for complete assessments. Our combination probe has a CPITN tip ideal for basic periodontal exams and a Williams for in-depth measurements. The probe markings are easily identifiable and in accordance with WHO guidelines.

**Front Face Mirror**

**Williams (10mm)/  
CPITN-C Probe 3.5/5.5/8.5/11.5**

150 Kits Per Box

**Kit: \$2.30**



# TWO PIECE DIAGNOSTIC KITS

BEST SELLER

‘Our front face mirror features ridges and grooves to improve grip stability, a zinc alloy tip for percussion testing and an anti-fog mirror surface.



**D2006**

## Diagnostic Kit

One of our best sellers, the combination perio probe/sickle is beneficial for both periodontal screening and caries detection.

**Front Face Mirror**

**Perio Probe 3/5/7/9/11/13/Sickle**

150 Kits Per Box

**Kit: \$2.00**



**BEST SELLER**

## TWO PIECE DIAGNOSTIC KITS



**D2009**

### Diagnostic Kit

Everything you need for a standard examination. New and improved, our front face mirror allows for high-quality illumination and retraction. Precision manufactured with stainless steel tips, our CPITN-C probe features easily identifiable markings in accordance with WHO guidelines.

**Front Face Mirror**

**Sickle/CPITN-C Probe 3.5/5.5/8.5/11.5**

150 Kits Per Box

**Kit: \$2.00**



# TWO PIECE DIAGNOSTIC KITS



**D2014**

## Diagnostic Kit

A variation of our ever popular 2 piece examination kit with a ball tip CPITN-C probe.

**Front Face Mirror**

**No. 8/CPITN-C Probe**

150 Kits Per Box

**Kit: \$2.00**



# THREE PIECE DIAGNOSTIC KITS

For a more detailed examination, check out our three-piece diagnostic kits. The additional tweezer or ruler gives you with that bit more insight into patient health.





# THREE PIECE DIAGNOSTIC KITS



**D3014**

## Diagnostic Kit

Featuring our new and improved front face mirror and examination tweezers as well as a CPITN-C ball tip probe

**Front Face Mirror**

**No. 8/CPITN-C Probe**

**Examination Tweezers**

100 Kits Per Box

**Kit: \$2.20**



**BEST SELLER**

## THREE PIECE DIAGNOSTIC KITS

**'Strong and sturdy with enhanced grip features on the handle and tips. Our premium London College Tweezers are an essential tool in most procedures.'**



**D3001**

### Premium Diagnostic Kit

Our D3001 is flying off the shelves for good reason. Featuring our London college tweezers, this three-piece combination assists in delivering premium quality examinations.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

London College Tweezers

100 Kits Per Box

**Kit: \$4.30**



## THREE PIECE DIAGNOSTIC KITS

**'When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.'**



**D3003**

### Orthodontic Kit

Our orthodontic exam kits will be your go-to for quick assessments. The small ruler size is ideal for both paediatric and adult consultations.

**Front Face Mirror**

**Perio Probe 3/5/7/9/11/13/Sickle**

**Orthodontic Ruler**

100 Kits Per Box

**Kit: \$2.30**



# THREE PIECE DIAGNOSTIC KITS



**D3004**

## Orthodontic Kit

Featuring an orthodontic ruler, our D3004 is the kit of choice for screening patients prior to referring them to the Orthodontist.

**Front Face Mirror**

**No. 8 Standard/CPITN-C Probe 3.5/5.5/8.5/11.5**

**Orthodontic Ruler**

100 Kits Per Box

**Kit: \$2.30**





# THREE PIECE DIAGNOSTIC KITS



**D3005**

## Premium Diagnostic Kit

In one simple step you will be ready for your patient. Arriving sterile and ready to use, patient screening has never been this simple.

Front Face Mirror

No. 8 Standard/CPITN-C Probe 3.5/5.5/8.5/11.5

London College Tweezers

100 Kits Per Box

**Kit: \$4.30**





**BEST SELLER**

## THREE PIECE DIAGNOSTIC KITS

**'Our economy tweezer features serrated, non-slip tips to ensure you have complete control during procedures.'**



**D3006**

### Diagnostic Kit

With a new and improved mirror and examination tweezer you can be confident in the quality of our D3006.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

Examination Tweezers

100 Kits Per Box

**Kit: \$2.20**



## THREE PIECE DIAGNOSTIC KITS



**'Looking for a ball tip probe?  
Our CPITN-C combinations have  
a 0.5mm ball tip and markings in  
accordance with WHO guidelines.'**

**D3010**

### Premium Diagnostic Kit

Our improved mirror design features a zinc alloy tip for percussion testing, making the D3010 ideal for indicating extraction or root treatment.

Front Face Mirror

Sickle/CPITN-C Probe 3.5/5.5/8.5/11.5

London College Tweezers

100 Kits Per Box

**Kit: \$4.30**





## **BEST SELLER**

**Say goodbye to the days of scrubbing mixing materials from your instruments with our S2002. This single-use mixing kit gives you spotless tools every time.**

# COMPLETE PERIODONTAL DIAGNOSTIC KIT

When a thorough assessment is required, our complete periodontal kit is for you.





**BEST SELLER**

## PERIODONTAL DIAGNOSTIC KITS



**'Our Waste to Energy service allows your clinic to reduce instrument waste going to landfill by 90%.'**

**SI550ST**

### Complete Diagnostic Kit

Convenient for mobile dentists and domiciliary care. With everything at your fingertips, our diagnostic pack is making patient examinations quick and easy.

25 Kits Per Box

**Kit: \$3.56**

- Front Face Mirror
- Perio Probe 3/5/7/9/11/13/Sickle
- Examination Tweezers
- 2 Face Masks with Ear Loops
- Patient Bib with Ties
- 2x Cotton Roll No. 2
- Tissue wipe
- Tray



# PACKAGING

Our instruments come conveniently packaged in sterile plastic sleeves, with clear viewing panels and concise labelling for easy identification.

The individual instruments and kit combinations are contained in front opening dispensers to provide swift access when stored in supply cupboards, on benchtops or shelves. These designs feature alternating images of our best selling instruments and kits, should you need some inspiration for your next purchase.

All levels of our product packaging come with the addition of **GS1 Data Matrix labels**. Switching to the standardised electronic communication system provides complete traceability of our instruments from the point of manufacture to the point of care.

We have worked hard to ensure that all our products and packaging are composed of recyclable materials.



# MATERIALS

## Instruments

Our instrument handles are Acrylonitrile Butadiene Styrene (ABS) plastic.

The instrument stems and tips are made from Stainless Steel, including our London College Tweezers. Our Examination Tweezers are Steel and Zinc coated.

The coating of our Front Face Mirror is Aluminium, the substrate is Polycarbonate (PC) and the Percussion Tip is Zinc Alloy.

Our Orthodontic Ruler is made from Polypropylene (PP).

In our Complete Periodontal Diagnostic Kits, the Tray is Polypropylene (PP), the Face Masks are a Non-Woven Fabric, the Patient Bib is Paper coated with Polyethylene film, the Tissue Wipe is Paper and the Cotton Rolls are... Cotton.

## Sleeves

Our sleeve design comprises of Polyethylene Terephthalate (PET) with a Tyvek® medical dialysing paper strip.

## Dispensers and Outer Cartons

The dispenser in which your instruments and kits are stored and the outer carton in which the products are shipped are both made up of 100% recycled paper and cardboard.

## Sharps Containers

The Terumo containers we supply are manufactured from high quality, puncture resistant Polypropylene (PP) and Polyethylene (PE).

All the above materials are recyclable.



# HOW TO BUY

## Online Ordering

When you are ready to place an order simply go to our website here, register and create an account and purchase the items you are interested in. Once an account is created you will then have your 'My Account' dashboard where you can view past orders, address details and account information. You can also choose to apply for a Trade Account if this is what you would prefer.

## Samples

We offer free samples to those who have not yet tried our instruments and would like to 'try before you buy'. If you would like to trial any instruments prior to purchasing please contact us.



# SHARPS CONTAINERS

As authorised distributors of Terumo sharps containers, we are now offering a clinical SMART Waste Recovery Program where we can dispose of your sharps waste.

With the aim of simplifying waste recovery forecasting for each clinic, the range of containers we have selected to supply hold similar capacities to our instrument boxes. An 11L container can dispose of around 1 dispenser, which is about 300 instruments. A 17L container fits around 2 dispensers, 600 instruments.

**11L Container with Roll Top Lid**  
including supply, collection & disposal

**17L Container with Roll Top Lid**  
including supply, collection & disposal

Easy to assemble, the containers are supplied nested to save space and with a separate lid that must be pressed down until all the locking clips engage. Once it's full it can be permanently locked and stored ready for collection. The accessories necessary to install the sharps containers in your clinic are also available.

Waste management service not yet available in all areas of Australia. Please contact our team or your state brand representative for more information.



**Single Use Dental Instruments**  
**1/24 Hoopers Road, Kunda Park**  
**QLD Australia 4556**  
**07 5445 6486**  
**[sales@sudinstruments.com](mailto:sales@sudinstruments.com)**  
**[www.singleusedentalinstruments.com](http://www.singleusedentalinstruments.com)**

