

**SINGLE USE**



**RESTORATIVE  
Catalogue**

# KEEPING IT SIMPLE



- > PRECISE
- > QUALITY
- > ERGONOMIC
- > STERILE
- > SINGLE USE

# CONTENTS

- 4 ABOUT US
- 5 OUR PRODUCTS
- 6 WASTE MANAGEMENT
- 8 TWO PIECE  
RESTORATIVE KITS
- 14 THREE PIECE  
RESTORATIVE KITS
- 18 FOUR PIECE  
RESTORATIVE KITS
- 24 FIVE PIECE  
RESTORATIVE KITS
- 28 SIX PIECE  
RESTORATIVE KITS
- 32 PACKAGING
- 33 MATERIALS
- 34 HOW TO BUY
- 35 SHARPS CONTAINERS

# ABOUT US

Our goal is to simplify dentistry as you know it. Minimising instrument preparation into one easy step, our single-use instruments will streamline the way you work. Quick, convenient and stress-free, they are improving dental clinics on a daily basis.

Pre-packed and sterile, they arrive ready for use. Removing the need for sterilisation allows you to focus on what's important: delivering high-quality patient care.

We know running an efficient clinic can be a juggling act. Our instruments will save your team time and take the stress out of a busy schedule.

Join us and take a convenient approach to dentistry.



# OUR PRODUCTS

By listening to customers and continuously working to improve our offering, our range has evolved to feature a wide variety of instruments and kits ideal for every procedure.

## Individual Instruments

Using individually wrapped instruments is a great way to transition your clinic towards implementing single-use solutions. Only open those you need for procedures and take minimising waste to a whole new level.

## Diagnostic Kits

Have you been entertaining the idea of single-use instruments? Equipping your clinic with our readily available diagnostic kits allows your nurse to spend less time setting up and more time chair side, focusing on the patient. Who would complain about that?

## Restorative Kits

Our restorative instrument kits will leave you staring at the autoclave wondering why you spend so much money maintaining it every year. Using our thoughtfully designed combinations alleviates the pressure and expense of reprocessing. Remember if you can't find a restorative combination that works for you, we'll make you one. That simple!



# WASTE MANAGEMENT

Our disposable instruments solve many problems for the dental industry. Reduced exposure to infection, minimised reprocessing and complete convenience... they really are the obvious choice.

Our mission is to ensure that our instruments become a single-use, zero waste solution. Through determination and collaboration with stakeholders who have a similar goal we have made significant progress in recent years.

The current waste management program is focused on a method that utilises a neutral waste recycling solution. The process minimises carbon emissions and instead repurposes the instruments into an energy source in the form of electricity, heat or fuel.

As we implement this proposal and develop a better understanding of our waste potential, we will be able to utilise this information to develop relationships with companies that want to recover and repurpose the raw material resources. This is our future vision and goal.






We are now offering a SMART Waste Recovery Program which enables clinics to dispose of their instrument waste into an allocated sharps container, have it collected, then taken to be processed within a dedicated waste facility.

To take advantage of our program all you need to do is register your interest via our website [here](#) or simply call our team on 07 5445 6486 and we can chat further to see if this program is suitable for your requirements. You can also find some more information [here](#).

The program includes the supply to clinics, collection when full and disposal of the waste.

See how the process works below.



\* Browse our sharps container range towards the back of this catalogue.

# TWO PIECE RESTORATIVE KITS

Thinking of phasing disposable instruments into your restorative procedures? Our thoughtfully designed two-piece kits can compliment many stages of a restoration.





## TWO PIECE RESTORATIVE KITS



**S2000**

### Restorative Kit

Suitable for any filling, this two-piece kit is a popular choice.

**Flat Plastic**

**Ball Burnisher**

150 Kits Per Box

**Kit: \$4.30**



# TWO PIECE RESTORATIVE KITS

BEST SELLER



**S2002**

## Restorative Kit

Say goodbye to the days of scrubbing mixing materials from your instruments with our S2002. This single-use mixing kit gives you spotless tools every time.

**Dycal Applicator**

**Spatula**

150 Kits Per Box

**Kit: \$3.50**



## TWO PIECE RESTORATIVE KITS



**S2003**

### Restorative Kit

This two-piece kit is popular for prepping and polishing a cavity restoration.

**Ball Burnisher**

**Small Spoon Excavator (1.2mm)**

150 Kits Per Box

**Kit: \$4.30**



## TWO PIECE RESTORATIVE KITS

‘Disposal of instruments immediately after use reduces exposure to infection and sharps injuries.’



**S2008**

### Restorative Kit

An ideal kit for placing and finishing Amalgam. Having a new plugger in every procedure means you never have to check for leftover residue after reprocessing. The small and large tips allow you to pack and contour different sized cavities.

**Amalgam Plugger Small/Large**

**Carver**

150 Kits Per Box

**Kit: \$4.30**



## TWO PIECE RESTORATIVE KITS



**S2012**

### Restorative Kit

Our S2012 includes a combination flat plastic/ball burnisher which means you can prep, place and polish with just two instruments.

Medium Spoon Excavator (1.8mm)

Flat Plastic/Ball Burnisher

150 Kits Per Box

**Kit: \$4.30**



# THREE PIECE RESTORATIVE KITS

Our three-piece kits are often our best selling combinations. Browse the selection below and see why they're a must have for so many of clinicians.





BEST SELLER

## THREE PIECE RESTORATIVE KITS



**S3000**

### Restorative Kit

One of our best sellers. With instruments for both placement and caries removal, this three-piece combination is perfect for most fillings.

Flat Plastic

Ball Burnisher

Medium Spoon Excavator (1.8mm)

100 Kits Per Box

**Kit: \$6.60**





# THREE PIECE RESTORATIVE KITS



‘The majority of our instruments are doubled-ended, giving you that much more bang for your buck.’

S3012

## Restorative Kit

From packing to perfection this kit has you covered. Your Amalgam restorations have never been so convenient.


Round Plugger (0.8mm/1.0mm)

Amalgam Plugger Small/Large

Half Hollenback

100 Kits Per Box

Kit: \$6.60





# FOUR PIECE RESTORATIVE KITS

Be ready for whatever walks through the door with our grab and go four-piece combinations. No fuss after use, they'll have you ready to walk out on time too!



## FOUR PIECE RESTORATIVE KITS



**\$4000**

### Restorative Kit

This four-piece combination has everything you need to keep things clean and dry before placing composite. The inclusion of our premium London college tweezers makes this kit a favourite for restorative procedures.

**Flat Plastic**

**Ball Burnisher**

**Medium Spoon Excavator (1.8mm)**

**London College Tweezers**

80 Kits Per Box

**Kit: \$8.10**



# FOUR PIECE RESTORATIVE KITS

BEST SELLER

'Not sure what to purchase?  
Our dispensers feature  
alternating images of our  
best selling instruments  
and kits, should you need  
some inspiration.'



S4001

## Restorative Kit

One of our best sellers this kit features our unique combination flat plastic/ball burnisher, giving you two instruments in one. A great value kit which is endorsed by dentists throughout Australia, it is a go-to for your restorative needs.

Flat Plastic/Ball Burnisher


Ladmore Plugger

Small Spoon Excavator (1.2mm)

Medium Spoon Excavator (1.8mm)

80 Kits Per Box

Kit: \$8.10





# FOUR PIECE RESTORATIVE KITS



**S4003**

## Restorative Kit

The S4003 just made prepping, placing and polishing a whole lot quicker.

Flat Plastic

Ball Burnisher

Small Spoon Excavator (1.2mm)

Medium Spoon Excavator (1.8mm)

80 Kits Per Box

**Kit: \$8.10**



# FOUR PIECE RESTORATIVE KITS



**S4004**

## Restorative Kit

With the inclusion of our small spoon excavator, the S4004 is perfect for paediatric restorations.

Flat Plastic

Ball Burnisher

Small Spoon Excavator (1.2mm)

London College Tweezers

80 Kits Per Box

**Kit: \$8.10**







**When a thorough assessment is required, our complete periodontal diagnostic kits are for you.**

# FIVE PIECE RESTORATIVE KITS

Whether you're removing cavities, placing filling materials or carving dental fillings, we have the five-piece instrument kit for you.



## FIVE PIECE RESTORATIVE KITS

'Our stems and tips are comprised of stainless steel, which provide a strong, sturdy and reliable instrument.'



**S5000**

### Restorative Kit

For a more complete Amalgam kit, look no further than our S5000.

Flat Plastic

Ball Burnisher


Carver

Amalgam Plugger Small/Large

Medium Spoon Excavator (1.8mm)

60 Kits Per Box

**Kit: \$9.90**



# FIVE PIECE RESTORATIVE KITS



**S5002**

## Restorative Kit

The various sized excavators in our S5002 will have you prepared for any cavity.

Flat Plastic

Ball Burnisher


Small Spoon Excavator (1.2mm)

Medium Spoon Excavator (1.8mm)

Large Spoon Excavator (2.3mm)

60 Kits Per Box

**Kit: \$9.90**





# FIVE PIECE RESTORATIVE KITS



**S5005**

## Restorative Kit

You'll be flying through your restorations with our five-piece restorative kit.

Flat Plastic

Ball Burnisher


Carver

Ladmore Plugger

Medium Spoon Excavator (1.8mm)

60 Kits Per Box

**Kit: \$9.90**



# SIX PIECE RESTORATIVE KITS

Last minute appointment? You can set your tray up quickly with our comprehensive six-piece kits.



# SIX PIECE RESTORATIVE KITS



**'When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.'**

**S6000**


## Restorative Kit

The more complex the procedure the more instruments you require. Simplify your preparation process into one simple step. Our S6000 is ready when you are.

|                                |                 |
|--------------------------------|-----------------|
| Half Hollenback                | Flat Plastic    |
| Small Spoon Excavator (1.2mm)  | Ball Burnisher  |
| Medium Spoon Excavator (1.8mm) | PFI49 Burnisher |

50 Kits Per Box

**Kit: \$11.50**





# SIX PIECE RESTORATIVE KITS

BEST SELLER



**S6001**


## Restorative Kit

For complex restorations, our S6001 has everything you need. Reducing instrument preparation into one easy step, this kit allows you to spend your time focused on the patient.

- |                                |                |
|--------------------------------|----------------|
| Ladmore Plugger                | Flat Plastic   |
| Amalgam Plugger Small/Large    | Ball Burnisher |
| Medium Spoon Excavator (1.8mm) | Carver         |

50 Kits Per Box

**Kit: \$11.50**





**An essential item in any dental practice, our front face mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability.**

# PACKAGING

Our instruments come conveniently packaged in sterile plastic sleeves, with clear viewing panels and concise labelling for easy identification.

The individual instruments and kit combinations are contained in front opening dispensers to provide swift access when stored in supply cupboards, on benchtops or shelves. These designs feature alternating images of our best selling instruments and kits, should you need some inspiration for your next purchase.

All levels of our product packaging come with the addition of [GS1 Data Matrix labels](#). Switching to the standardised electronic communication system provides complete traceability of our instruments from the point of manufacture to the point of care.

We have worked hard to ensure that all our products and packaging are composed of recyclable materials.



# MATERIALS

## Instruments

Our instrument handles are Acrylonitrile Butadiene Styrene (ABS) plastic.

The instrument stems and tips are made from Stainless Steel, including our London College Tweezers. Our Examination Tweezers are Steel and Zinc coated.

The coating of our Front Face Mirror is Aluminium, the substrate is Polycarbonate (PC) and the Percussion Tip is Zinc Alloy.

Our Orthodontic Ruler is made from Polypropylene (PP).

In our Complete Periodontal Diagnostic Kits, the Tray is Polypropylene (PP), the Face Masks are a Non-Woven Fabric, the Patient Bib is Paper coated with Polyethylene film, the Tissue Wipe is Paper and the Cotton Rolls are... Cotton.

## Sleeves

Our sleeve design comprises of Polyethylene Terephthalate (PET) with a Tyvek® medical dialysing paper strip.

## Dispensers and Outer Cartons

The dispenser in which your instruments and kits are stored and the outer carton in which the products are shipped are both made up of 100% recycled paper and cardboard.

## Sharps Containers

The Terumo containers we supply are manufactured from high quality, puncture resistant Polypropylene (PP) and Polyethylene (PE).

All the above materials are recyclable.





# HOW TO BUY

## Online Ordering

When you are ready to place an order simply go to our website here, register and create an account and purchase the items you are interested in. Once an account is created you will then have your 'My Account' dashboard where you can view past orders, address details and account information. You can also choose to apply for a Trade Account if this is what you would prefer.

## Samples

We offer free samples to those who have not yet tried our instruments and would like to 'try before you buy'. If you would like to trial any instruments prior to purchasing please contact us.



# SHARPS CONTAINERS

As authorised distributors of Terumo sharps containers, we are now offering a clinical SMART Waste Recovery Program where we can dispose of your sharps waste.

With the aim of simplifying waste recovery forecasting for each clinic, the range of containers we have selected to supply hold similar capacities to our instrument boxes. An 11L container can dispose of around 1 dispenser, which is about 300 instruments. A 17L container fits around 2 dispensers, 600 instruments.

**11L Container with Roll Top Lid**  
including supply, collection & disposal

**17L Container with Roll Top Lid**  
including supply, collection & disposal

Easy to assemble, the containers are supplied nested to save space and with a separate lid that must be pressed down until all the locking clips engage. Once it's full it can be permanently locked and stored ready for collection. The accessories necessary to install the sharps containers in your clinic are also available.

Waste management service not yet available in all areas of Australia. Please contact our team or your state brand representative for more information.



**Single Use Dental Instruments**  
**1/24 Hoopers Road, Kunda Park**  
**QLD Australia 4556**  
**07 5445 6486**

**[sales@sudinstruments.com](mailto:sales@sudinstruments.com)**

**[www.singleusedentalinstruments.com](http://www.singleusedentalinstruments.com)**

