

SINGLE USE

INDIVIDUAL INSTRUMENT Catalogue



KEEPING IT SIMPLE

- > PRECISE
- > QUALITY
- > ERGONOMIC
- > STERILE
- > SINGLE USE



CONTENTS

- Click the links to the pages you want to view.

- [4 ABOUT US](#)
- [5 OUR PRODUCTS](#)
- [6 WASTE MANAGEMENT](#)
- [8 BURNISHERS](#)
- [10 CARVERS](#)
- [12 CURETTES](#)
- [15 EXCAVATORS](#)
- [17 FLAT PLASTICS](#)
- [19 MIRRORS](#)
- [20 PLUGGERS](#)
- [22 PROBES](#)
- [26 SCALERS](#)
- [28 TWEEZERS](#)
- [30 OTHER INSTRUMENTS](#)
- [32 PACKAGING](#)
- [33 MATERIALS](#)
- [34 HOW TO BUY](#)
- [35 SHARPS CONTAINERS](#)

ABOUT US

Our goal is to simplify dentistry as you know it. Minimising instrument preparation into one easy step, our single-use instruments will streamline the way you work. Quick, convenient and stress-free, they are improving dental clinics on a daily basis.

Pre-packed and sterile, they arrive ready for use. Removing the need for sterilisation allows you to focus on what's important: delivering high-quality patient care.

We know running an efficient clinic can be a juggling act. Our instruments will save your team time and take the stress out of a busy schedule.

Join us and take a convenient approach to dentistry.



OUR PRODUCTS

By listening to customers and continuously working to improve our offering, our range has evolved to feature a wide variety of instruments and kits ideal for every procedure.

Individual Instruments

Using individually wrapped instruments is a great way to transition your clinic towards implementing single-use solutions. Only open those you need for procedures and take minimising waste to a whole new level.

Diagnostic Kits

Have you been entertaining the idea of single-use instruments? Equipping your clinic with our readily available diagnostic kits allows your nurse to spend less time setting up and more time chair side, focusing on the patient. Who would complain about that?

Find out more in our [diagnostic range](#).

Restorative Kits

Our restorative instrument kits will leave you staring at the autoclave wondering why you spend so much money maintaining it every year. Using our thoughtfully designed combinations alleviates the pressure and expense of reprocessing. Remember if you can't find a restorative combination that works for you, we'll make you one. That simple!

View these kits in our [restorative range](#).



WASTE MANAGEMENT

Our disposable instruments solve many problems for the dental industry. Reduced exposure to infection, minimised reprocessing and complete convenience... they really are the obvious choice.

Our mission is to ensure that our instruments become a single-use, zero waste solution. Through sheer determination and collaboration with stakeholders who have a similar goal, we are much closer to making this a reality.

By utilising Waste to Energy technology, waste going to landfill is reduced by 90%. The process minimises carbon emissions and instead repurposes the instruments into an energy source in the form of electricity, heat or fuel.


We are now offering a waste collection service which enables clinics to dispose of their instrument waste into an allocated sharps container, have it collected, then taken to be processed within a dedicated waste facility.



To take advantage of our waste management service all you need to do is find the sharps container size that suits your clinic's needs and add it to your next instrument order. Once it reaches capacity, simply call our team on 07 5445 6486 to organise collection.

The price of the containers includes the cost of supply to clinics, collection when full and disposal of the waste. We do not believe in using better waste management as a tactic to increase profit. Our price is set to simply cover the cost of this improved process.

See how the process works below.



* Browse our sharps container range towards the back of this catalogue.

INDIVIDUAL INSTRUMENT RANGE

BURNISHERS

Our double-ended burnishers give you the benefit of two polishing tips in one. The different sized ends allow you to effectively adapt filling material with just one instrument.



BURNISHERS

- Click on the instrument image to view in more detail on our website.



BB49

PFI49 Burnisher

The PFI49 burnisher is ideal for placing and contouring restorations. Having both a rounded and pointed tip makes this instrument useful for placing composite in small increments prior to curing.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER



BB56

Ball Burnisher

The perfect blending tool. Our ball burnisher has both small and large tip ends to effectively condense and shape the composite against the cavity wall surfaces.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER



CB59

Flat Plastic/ Ball Burnisher


The combined flat plastic and ball burnisher is brilliant for placing and blending composite material.

20 Instruments Per Box

Instrument: \$2.25

CARVERS

You'll never feel the frustration of shaping a restoration with an instrument that's just not cutting it again. Our single-use carvers are guaranteed sharp and effective for all your contouring needs.



'When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.'

CARVERS

CV55

Carver

Being single-use, our carver tips are always sharp and perfectly pointed.

20 Instruments Per Box

Instrument: \$2.25

CV62

Wards Carver No.2

Sterile, sharp and ready for use, our CV62 is about to carve convenience into your restorations.

20 Instruments Per Box

Instrument: \$2.25

CV66

Half Hollenback

The CV66 is a versatile instrument which is very useful for placing, carving and contouring amalgam.

20 Instruments Per Box

Instrument: \$2.25

CB60

Flat Plastic/Carver

Our CB60 combines a flat plastic and carver into one instrument allowing you to conveniently prep, shape and contour filling material.


20 Instruments Per Box

Instrument: \$2.25



CURETTES

Get confused by curette tips? So do we. The thick brightly coloured handles on our disposable curettes not only feel great in your hand, they're easy to identify.



CURETTES



***GC11**

Gracey Curette 11/12

Using disposable curettes gives you the optimum working end in each appointment. Our GC11 is a great choice for finishing subgingival calculus removal.

20 Instruments Per Box

Instrument: \$2.50



***GC13**

Gracey Curette 13/14

Adopting the use of our single-use GC13 means the lower cutting edge will be prime for scaling, every time. The smooth edges on our curettes are perfect for removing subgingival calculus and getting into pockets without harming the patient.

20 Instruments Per Box

Instrument: \$2.50



***GC17**

Gracey Curette 17/18

Our disposable GC17 will always have the perfect angle for reaching distal surfaces of posterior teeth. The thick handle design provides you with comfort and accentuated grip throughout procedures.

20 Instruments Per Box

Instrument: \$2.50

***Limited availability.**

CURETTES

***UC02**

Universal Curette 2R/2L

Never have a dull cutting edge again when you make the switch to our single-use UC02. Fresh out the pack, they provide you with reliable quality in your scaling procedures.

20 Instruments Per Box

Instrument: \$2.50



'Our Waste to Energy service allows your clinic to reduce instrument waste going to landfill by 90%.'

***UC04**

Universal Curette 4R/4L

Sterile and ready for use, our UC04 is the convenient choice for incisors, canines and premolars.

20 Instruments Per Box

Instrument: \$2.50

***Limited availability.**

EXCAVATORS

No matter the size of the cavity we have the spoon for the task. Have more control of your manual excavation with sizes ranging from extra small to large.



EXCAVATORS

XC66

Extra Small Spoon Excavator (0.8mm)

The go-to instrument for clearing out the smallest micro-dental preparations.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER

XC67

Small Spoon Excavator (1.2mm)

One of our best sellers, our small spoon excavator is the top choice when you want to have extra control over caries removal in deep cavities.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER



XC58

Medium Spoon Excavator (1.8mm)

Gone are the days of dull and blunt instruments. Precision made, our single-use XC58 will have sharp edges each and every time.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER

XC68

Large Spoon Excavator (2.3mm)

The large curved surface on our XC68 excavator is ideal for removing broad-scale decay.

20 Instruments Per Box

Instrument: \$2.25

FLAT PLASTICS

Small, regular or combination flat plastic?
Whatever your preference we have the one for you.



FLAT PLASTICS



FP48

Small Flat Plastic

About half the size of our original flat plastic, the tips on our FP48 are ideal for small and narrow cavities.

20 Instruments Per Box

Instrument: \$2.25



FP54

Flat Plastic

Our double-ended flat plastic features long, flat tips that are helpful for carrying and shaping composite materials.

20 Instruments Per Box

Instrument: \$2.25



CB59

Flat Plastic/ Ball Burnisher

The combined flat plastic and ball burnisher is brilliant for placing and blending composite material.

20 Instruments Per Box

Instrument: \$2.25



CB60

Flat Plastic/Carver

Our CB60 combines a flat plastic and carver into one instrument allowing you to conveniently prep, shape and contour filling material.

20 Instruments Per Box

Instrument: \$2.25

BEST SELLER

MIRROR

An essential item in any dental practice, our front face mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability. With the addition of a zinc alloy tip for percussion testing and an anti-fog mirror surface...

What more could you ask for?



'All our instruments and packaging are made from recyclable materials.'

FF53

Front Face Mirror

20 Instruments Per Box

Instrument: \$0.75

BEST SELLER

PLUGGERS

Our single-use pluggers can bring a new level of convenience to your endodontic and restorative treatments.



PLUGGERS

PL57

Ladmore Plugger

Our ladmore plugger features both medium and large tips which are ideal for placing and adapting composite.

20 Instruments Per Box

Instrument: \$2.25

PL61

Amalgam Plugger

Small/Large

The single-use amalgam plugger has a large and a small plugging surface. The combination is useful for both broad and deep amalgam condensing to ensure a strong filling.

20 Instruments Per Box

Instrument: \$2.25

PL62

Round Plugger (0.8mm/1.0mm)

The round plugger has two different sized ends for compaction from the base to the surface of the tooth. Perfectly smooth tips for every use allowing for a flawless finish.

20 Instruments Per Box

Instrument: \$2.25

PL63

Machtou Plugger (0.5mm/0.6mm)

An endodontic plugger with varied tip sizes for condensing material into different sized root canals.

20 Instruments Per Box

Instrument: \$2.25

PL64

Machtou Plugger (0.8mm/1.0mm)

Our manual endodontic pluggers are long and flexible for manoeuvrability within narrow and curved canals. The tips on our PL64 cater for larger canals.

20 Instruments Per Box

Instrument: \$2.25



PROBES

We stock a variety of disposable probe combination instruments which can help make the assessment and diagnosis of patient health that much quicker and easier in your clinic. Browse our range and see which one suits your style.



**‘Looking for a ball tip probe?
Our CPITN-C combinations have
a 0.5mm ball tip and markings in
accordance with WHO guidelines.’**

PROBES



PB84

No. 8 Standard/ CPITN-C Probe

Ensuring that our product range caters for everyone our PB84 is popular amongst dentists trained within the UK.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.30



PB85

Sickle/CPITN-C Probe

Our PB85 combines an explorer and CPITN-C probe into one instrument. Featuring precision edges and easily identifiable measurements, this probe is a convenient addition to diagnostic assessments.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.30



PB89

Perio Probe/Sickle

A customer favourite, our perio probe/sickle has easily identifiable markings to assist in delivering accurate assessments. This combination instrument can be purchased individually or within our top selling three-piece diagnostic kits: D3001 and D3006.

Perio Probe: 3/5/7/9/11/13mm.

20 Instruments Per Box

Instrument: \$1.30

BEST SELLER

PROBES



PB90

Williams (10mm)/ CPITN-C Probe

The PB90 features a Williams tip with the 4mm and 6mm markings omitted for easy identification and a CPITN-C probe with additional 8.5–11.5mm markings.

Williams: 1/2/3/5/7/8/9/10mm.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.60



PB91

UNC-15/CPITN-C Probe

Combining the popular ball tip/WHO probe with our UNC-15, this probe is an excellent choice when measuring pocket depth.

UNC-15: 1/2/3/4/5/6/7/8/9/10/11/12/13/14/15.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.70



PB92

No. 12/CPITN-C Probe

Our sterile combination probes make patient examinations quick and convenient. The ball tip maximises patient safety when measuring pockets.

No.12: 3/6/9/12mm.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.30



PB93

No. 6 Probe

Our single-ended No. 6 Probe is useful for many purposes. Whether you're inspecting filling margins or checking cavities prior to filling, this probe is in prime condition and ready to go.

20 Instruments Per Box

Instrument: \$1.00



PB94

No. 6/CPITN-C Probe

Our No. 6 and ball tip combination instrument allows you to probe the entire mouth without causing patient discomfort.

CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip.

20 Instruments Per Box

Instrument: \$1.30



PB95

No. 6/Williams Probe (10mm)

Probe, explore and measure with one instrument when using our combination

No. 6/Williams Probe.

Williams: 1/2/3/5/7/8/9/10mm.

20 Instruments Per Box

Instrument: \$1.60

SCALERS

Sharpening will be a thing of the past once you experience our single-use scalers. They're perfectly precise, every time.



'Minimising instrument preparation into one easy step lets you spend less time setting up and more time focusing on the patient.'

SCALERS

SC02

Sickle Scaler 204S

A posterior sickle scaler designed to remove supragingival interproximal calculus. Newly improved, our SC02 features a thick handle with ridges and grooves for improved comfort and grip.

20 Instruments Per Box

Instrument: \$2.50

SC03

Hygienist Scaler H6/H7

A favourite amongst hygienists this scaler is designed for interproximal cleaning between anterior teeth and premolars. Featuring a stainless-steel stem and ergonomic handle to ensure stability and control during procedures.

20 Instruments Per Box

Instrument: \$2.50

BEST SELLER



TWEEZERS

Whether you're shopping on a budget or looking for a premium option, we have the tweezer for you. Our Examination Tweezers are helpful in any procedure and our London College Tweezers are ideal when additional strength and grip is required.



TWEEZERS

TW82

London College Tweezers

Strong and sturdy with enhanced grip features on the handle and tips. Our premium London College Tweezers are an essential tool in most procedures.

20 Instruments Per Box

Instrument: \$2.50

BEST SELLER



TW81

Examination Tweezers

Our economy tweezer is a popular choice for both diagnostic and restorative consultations. Featuring serrated, non-slip tips to ensure you have complete control during procedures.

20 Instruments Per Box

Instrument: \$0.80

OTHER INSTRUMENTS



'Our brand representatives are available to tailor our product offering to suit your clinic's needs. Contact our team and see what we can do for you.'

OTHER INSTRUMENTS

DA69

Dycal Applicator

Our single-ended dycal applicator is strong and slender to ensure accurate placement of lining materials.

20 Instruments Per Box

Instrument: \$1.00

BEST SELLER



SP65

Spatula

The single-use spatula is slim, sleek and refined. Perfectly clean every time, these spatulas remove the reprocessing hassle after mixing and spreading sticky dental materials.

20 Instruments Per Box

Instrument: \$2.50

BEST SELLER



PACKAGING

Our instruments come conveniently packaged in sterile plastic sleeves, with clear viewing panels and concise labelling for easy identification.

The individual instruments are contained in a small and convenient box that provides swift access when stored in supply cupboards, on benchtops or shelves. This design features a range of our instruments, should you need some inspiration for your next purchase.

All levels of our product packaging come with the addition of [GS1 Data Matrix labels](#). Switching to the standardised electronic communication system provides complete traceability of our instruments from the point of manufacture to the point of care.

We have worked hard to ensure that all our products and packaging are composed of recyclable materials.



MATERIALS

Instruments

Our instrument handles are Acrylonitrile Butadiene Styrene (ABS) plastic.

The instrument stems and tips are made from Stainless Steel, including our London College Tweezers. Our Examination Tweezers are Steel and Zinc coated.

The coating of our Front Face Mirror is Aluminium, the substrate is Polycarbonate (PC) and the Percussion Tip is Zinc Alloy.

Our Orthodontic Ruler is made from Polypropylene (PP).

In our Complete Periodontal Diagnostic Kits, the Tray is Polypropylene (PP), the Face Masks are a Non-Woven Fabric, the Patient Bib is Paper coated with Polyethylene film, the Tissue Wipe is Paper and the Cotton Rolls are... Cotton.

Sleeves

Our sleeve design comprises of Polyethylene Terephthalate (PET) with a Tyvek® medical dialysing paper strip.

Boxes and Outer Cartons

The boxes in which your instruments and kits are stored and the outer carton in which the products are shipped are both made up of 100% recycled paper and cardboard.

Sharps Containers

The Terumo containers we supply are manufactured from high quality, puncture resistant Polypropylene (PP) and Polyethylene (PE).

All the above materials are recyclable.



HOW TO BUY

Online Quoting System

If you would like to obtain an accurate quote inclusive of GST and the delivery charge to your area, our newly launched website features an instant quoting system. All you have to do is visit our [product page](#), add your preferred items to your quote and fill in the required information to receive it instantly.

To make this quicker in future, you can [create an account](#) so you can easily review previous quotes and request new ones without re-entering your details.

Purchasing

When you are ready to place an order, simply contact our team via phone on 07 5445 6486 or [email](#) to purchase items you have been quoted for previously, or to request new ones. Our online quotes are valid for 30 days and can be emailed to our team directly for efficient, accurate order placement.

Samples

If you would like to trial any instruments prior to purchasing, free sterile samples can be requested [here](#).



SHARPS CONTAINERS

As authorised distributors of Terumo sharps containers, we are now offering a clinical waste management service where we can recycle your sharps waste into energy. When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.

*Waste management service not yet available in all areas of Australia. Please contact our team for more information.

With the aim of simplifying waste collection forecasting for each clinic, the range of containers we have selected to supply hold similar capacities to our instrument boxes. A single 11L waste container allows for disposal of approximately 300 instruments. The 17L waste container allows for disposal of approximately 600 instruments.

11L Container with Roll Top Lid
including supply, collection & disposal - \$28.00

17L Container with Roll Top Lid
including supply, collection & disposal - \$35.00

The accessories necessary to install the sharps containers in your clinic are also available. Please [contact](#) our team or your state brand representative for more information.



Single Use Dental Instruments
1/24 Hoopers Road, Kunda Park
QLD Australia 4556
07 5445 6486

sales@sudinstruments.com

www.singleusedentalinstruments.com

