

SINGLE USE



D3006

Featuring our new and improved Front Face Mirror, Perio Probe/Sickle and Examination Tweezers, our high-quality D3006 is a customer favourite.

Single Use Dental Instruments Pty Ltd

1/24 Hoopers Road, Kunda Park
QLD 4556 Australia
sales@sudinstruments.com
07 5445 6486

Please see website link for kit details [here](#).

ERGONOMIC • PRECISE • QUALITY • STERILE • SINGLE USE

FF53

Front Face Mirror


Description:

An essential item in any dental practice our Front Face Mirror is new and improved. With a focus on comfort, our ergonomic handle features additional ridges and grooves to improve grip stability. With the addition of a zinc alloy tip for percussion testing and an anti-fog mirror surface.

Features:

- Zinc alloy tip for percussion testing
- Ergonomic handle with ridges and grooves for improved grip
- Latex-free
- Arrives sterile and ready for use
- GS1 Data Matrix
- CE and TGA Certified

Materials: Zinc Alloy, Acrylonitrile butadiene styrene (ABS), Polycarbonate (PC), Aluminiumbutadiene styrene (ABS)



PB89

Perio Probe/Sickle


Description:

A customer favourite, our perio probe/sickle has easily identifiable markings to assist in delivering accurate assessments. This combination instrument can be purchased individually or within our top selling diagnostic kits; D3001 and D3006.

Features:

- Combination instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Calibrations; Perio Probe: 3/5/7/9/11/13mm
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



TW81

Examination Tweezers


Description:

Our economy Tweezer is a popular choice for both diagnostic and restorative consultations. Featuring serrated, non-slip tips to ensure you have complete control during procedures.

Features:

- Ergonomic handle
- Serrated non-slip tips
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Steel plate, zinc coating



Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

Quality assurance

All our instrument manufacturing has CE, TGA and ISO certification to ensure quality compliance with industry standards.

Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

