

SINGLE USE



S3012

From packing to perfection,
this kit has you covered.
Your Amalgam restorations
have never been so
convenient with our 1/2
Hollenback, Amalgam
Plugger (Large/Small)
and Round Plugger
(0.8/1.0mm).

Please see website link for kit details [here](#).

Single Use Dental Instruments Pty Ltd

1/24 Hoopers Road, Kunda Park
QLD 4556 Australia
sales@sudinstruments.com
07 5445 6486

ERGONOMIC • PRECISE • QUALITY • STERILE • SINGLE USE

CV66

1/2 Hollenback


Description:

Our 1/2 Hollenback is a versatile instrument used to place, carve and contour amalgam.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



PL61

Amalgam Plugger (Small/Large)


Description:

The single-use Amalgam Plugger has both a large and a small plugging surface. The combination is useful for both broad and deep amalgam condensing to ensure a high-quality filling.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 2.5mm/1.5mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



PL62

Round Plugger 0.8mm/1.0mm


Description:

The Round Plugger has two different sized ends for compaction from the base to the surface of the tooth. Perfectly smooth tips for every use allowing for a flawless finish.

Features:

- Double-ended instrument
- Precision manufactured stainless steel tips
- Tip diameters; 0.8mm/1.0mm
- Ergonomic handle with ridges and grooves for improved grip
- Stainless steel stem for added stability
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE and TGA Certified

Materials: Stainless steel, Acrylonitrile butadiene styrene (ABS)



Our Products

Our unique, single-use instruments were designed with simplicity in mind. Featuring a broad variety of diagnostic and restorative instruments, our goal is to remove daily stress for clinicians and improve patient safety. Not only convenient, we have created a product range that is focused on quality and sustainability. Our single use instruments arrive sterile and ready for immediate use. The ergonomic handle is made from ABS plastic and features ridges and grooves for an improved grip. The stainless-steel stem provides additional weight and durability to ensure it is fit for purpose. Precision manufactured stainless steel tips replicate the design and function of re-usable instruments.

Packaging and materials

We are committed to reducing the environmental impact of our products. As such our products and packaging are manufactured with recyclable materials. Additionally, we source recycled cardboard to create our benchtop dispensers and shipping containers.

Quality assurance

All our instrument manufacturing has CE, TGA and ISO certification to ensure quality compliance with industry standards.

Tyvek® medical dialysing strip

A Tyvek® medical dialysing strip is incorporated into all our product sleeves. Comprising of flashspun and bonded high-density polyethylene (HDPE) filaments, the strip provides a microbial barrier and exceptional tear and puncture resistance to maintain product integrity during transit.

