

SINGLE USE



DIAGNOSTIC Catalogue

KEEPING IT SIMPLE

- > PRECISE
- > QUALITY
- > ERGONOMIC
- > STERILE
- > SINGLE USE



CONTENTS

- Click the links to the pages you want to view.

- 4 [ABOUT US](#)
- 5 [OUR PRODUCTS](#)
- 6 [WASTE MANAGEMENT](#)
- 8 [TWO PIECE
DIAGNOSTIC KITS](#)
- 12 [THREE PIECE
DIAGNOSTIC KITS](#)
- 20 [COMPLETE PERIODONTAL
DIAGNOSTIC KITS](#)
- 24 [PACKAGING](#)
- 25 [MATERIALS](#)
- 26 [HOW TO BUY](#)
- 27 [SHARPS CONTAINERS](#)

ABOUT US

Our goal is to simplify dentistry as you know it. Minimising instrument preparation into one easy step, our single-use instruments will streamline the way you work. Quick, convenient and stress-free, they are improving dental clinics on a daily basis.

Pre-packed and sterile, they arrive ready for use. Removing the need for sterilisation allows you to focus on what's important: delivering high-quality patient care.

We know running an efficient clinic can be a juggling act. Our instruments will save your team time and take the stress out of a busy schedule.

Join us and take a convenient approach to dentistry.



OUR PRODUCTS

By listening to customers and continuously working to improve our offering, our range has evolved to feature a wide variety of instruments and kits ideal for every procedure.

Individual Instruments

Using individually wrapped instruments is a great way to transition your clinic towards implementing single-use solutions. Only open those you need for procedures and take minimising waste to a whole new level.

Find out more by viewing our [individual range](#).

Diagnostic Kits

Have you been entertaining the idea of single-use instruments? Equipping your clinic with our readily available diagnostic kits allows your nurse to spend less time setting up and more time chair side, focusing on the patient. Who would complain about that?

Restorative Kits

Our restorative instrument kits will leave you staring at the autoclave wondering why you spend so much money maintaining it every year. Using our thoughtfully designed combinations alleviates the pressure and expense of reprocessing. Remember if you can't find a restorative combination that works for you, we'll make you one. That simple!

View these kits in our [restorative range](#).



WASTE MANAGEMENT

Our disposable instruments solve many problems for the dental industry. Reduced exposure to infection, minimised reprocessing and complete convenience... they really are the obvious choice.

Our mission is to ensure that our instruments become a single-use, zero waste solution. Through sheer determination and collaboration with stakeholders who have a similar goal, we are much closer to making this a reality.

By utilising Waste to Energy technology, waste going to landfill is reduced by 90%. The process minimises carbon emissions and instead repurposes the instruments into an energy source in the form of electricity, heat or fuel.


We are now offering a waste collection service which enables clinics to dispose of their instrument waste into an allocated sharps container, have it collected, then taken to be processed within a dedicated waste facility.



To take advantage of our waste management service all you need to do is find the sharps container size that suits your clinic's needs and add it to your next instrument order. Once it reaches capacity, simply call our team on 07 5445 6486 to organise collection.

The price of the containers includes the cost of supply to clinics, collection when full and disposal of the waste. We do not believe in using better waste management as a tactic to increase profit. Our price is set to simply cover the cost of this improved process.

See how the process works below.



* Browse our sharps container range towards the back of this catalogue.

TWO PIECE DIAGNOSTIC KITS

Our two-piece diagnostic kits are ideal for an initial check of the mouth. The anti-fog mirror delivers clear vision and our combined probes give you two instruments in one. Browse our selection to see which mirror and probe combination suits you best.



TWO PIECE DIAGNOSTIC KITS



- Click the 'view more images' button to see each kit in more detail on our website.



D2004

Diagnostic Kit

Great for complete assessments. Our combination probe has a CPITN tip ideal for basic periodontal exams and a Williams for in-depth measurements. The probe markings are easily identifiable and in accordance with WHO guidelines.

Front Face Mirror

**Williams (10mm)/
CPITN-C Probe 3.5/5.5/8.5/11.5**

150 Kits Per Box

Kit: \$2.30

VIEW MORE IMAGES

TWO PIECE DIAGNOSTIC KITS

BEST SELLER

Our front face mirror features ridges and grooves to improve grip stability, a zinc alloy tip for percussion testing and an anti-fog mirror surface.



D2006

Diagnostic Kit

One of our best sellers, the combination perio probe/sickle is beneficial for both periodontal screening and caries detection.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

150 Kits Per Box

Kit: \$2.00



[VIEW MORE IMAGES](#)

BEST SELLER

TWO PIECE DIAGNOSTIC KITS



D2009

Diagnostic Kit

Everything you need for a standard examination. New and improved, our front face mirror allows for high-quality illumination and retraction. Precision manufactured with stainless steel tips, our CPITN-C probe features easily identifiable markings in accordance with WHO guidelines.

Front Face Mirror

Sickle/CPITN-C Probe 3.5/5.5/8.5/11.5

150 Kits Per Box

Kit: \$2.00



[VIEW MORE IMAGES](#)

THREE PIECE DIAGNOSTIC KITS

For a more detailed examination, check out our three-piece diagnostic kits. The additional tweezer or ruler gives you with that bit more insight into patient health.



BEST SELLER

THREE PIECE DIAGNOSTIC KITS

'Strong and sturdy with enhanced grip features on the handle and tips. Our premium London College Tweezers are an essential tool in most procedures.'



D3001

Premium Diagnostic Kit

Our D3001 is flying off the shelves for good reason. Featuring our London college tweezers, this three-piece combination assists in delivering premium quality examinations.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

London College Tweezers

100 Kits Per Box

Kit: \$4.30



VIEW MORE IMAGES

THREE PIECE DIAGNOSTIC KITS

'When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.'



D3003

Orthodontic Kit

Our orthodontic exam kits will be your go-to for quick assessments. The small ruler size is ideal for both paediatric and adult consultations.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

Orthodontic Ruler

100 Kits Per Box

Kit: \$2.30



[VIEW MORE IMAGES](#)

THREE PIECE DIAGNOSTIC KITS



D3004

Orthodontic Kit

Featuring an orthodontic ruler, our D3004 is the kit of choice for screening patients prior to referring them to the Orthodontist.

Front Face Mirror

No. 8 Standard/CPITN-C Probe 3.5/5.5/8.5/11.5

Orthodontic Ruler

100 Kits Per Box

Kit: \$2.30



[VIEW MORE IMAGES](#)

THREE PIECE DIAGNOSTIC KITS



D3005

Premium Diagnostic Kit

In one simple step you will be ready for your patient. Arriving sterile and ready to use, patient screening has never been this simple.

Front Face Mirror

No. 8 Standard/CPITN-C Probe 3.5/5.5/8.5/11.5

London College Tweezers

100 Kits Per Box

Kit: \$4.30



[VIEW MORE IMAGES](#)

BEST SELLER

THREE PIECE DIAGNOSTIC KITS

‘Our economy tweezer features serrated, non-slip tips to ensure you have complete control during procedures.’



D3006

Diagnostic Kit

With a new and improved mirror and examination tweezer you can be confident in the quality of our D3006.

Front Face Mirror

Perio Probe 3/5/7/9/11/13/Sickle

Examination Tweezers

100 Kits Per Box

Kit: \$2.20

[VIEW MORE IMAGES](#)

THREE PIECE DIAGNOSTIC KITS



**'Looking for a ball tip probe?
Our CPITN-C combinations have
a 0.5mm ball tip and markings in
accordance with WHO guidelines.'**

D3010

Premium Diagnostic Kit

Our improved mirror design features a zinc alloy tip for percussion testing, making the D3010 ideal for indicating extraction or root treatment.

Front Face Mirror

Sickle/CPITN-C Probe 3.5/5.5/8.5/11.5

London College Tweezers

100 Kits Per Box

Kit: \$4.30



[VIEW MORE IMAGES](#)



BEST SELLER

Say goodbye to the days of scrubbing mixing materials from your instruments with our [S2002](#). This single-use mixing kit gives you spotless tools every time.

[BROWSE OUR RESTORATIVE RANGE HERE](#)

COMPLETE PERIODONTAL DIAGNOSTIC KITS

When a thorough assessment is required, our complete periodontal kits are for you. From mirrors to masks and the option between two different probes, they have everything you need for patient examinations.



BEST SELLER

PERIODONTAL DIAGNOSTIC KITS



'Our Waste to Energy service allows your clinic to reduce instrument waste going to landfill by 90%.'

- Front Face Mirror
- Perio Probe 3/5/7/9/11/13/Sickle
- Examination Tweezers
- 2 Face Masks with Ear Loops
- Patient Bib with Ties
- 2x Cotton Roll No. 2
- Tissue wipe
- Tray

SI550ST

Complete Diagnostic Kit

Convenient for mobile dentists and domiciliary care. With everything at your fingertips, our diagnostic pack is making patient examinations quick and easy.

25 Kits Per Box

Kit: \$2.70

VIEW MORE IMAGES

PERIODONTAL DIAGNOSTIC KITS



SI950ST

Complete Diagnostic Kit

Exams on the go?

Our comprehensive examination kit is easy to transport and has everything you need for patient screening.

25 Kits Per Box

Kit: \$2.80

Front Face Mirror
No. 8 Standard/CPITN-C Probe 3.5/5.5/8.5/11.5
Examination Tweezers
2 Face Masks with Ear Loops
Patient Bib with Ties
2x Cotton Roll No. 2
Tissue wipe
Tray

[VIEW MORE IMAGES](#)



BEST SELLER

A customer favourite, our PB89 perio probe/sickle has easily identifiable markings to assist in delivering accurate assessments. This can be purchased individually or within our top selling three-piece diagnostic kits: [D3001](#) and [D3006](#).

[BROWSE OUR INDIVIDUAL RANGE HERE](#)

PACKAGING

Our instruments come conveniently packaged in sterile plastic sleeves, with clear viewing panels and concise labelling for easy identification.

The individual instruments and kit combinations are contained in front opening dispensers to provide swift access when stored in supply cupboards, on benchtops or shelves. These designs feature alternating images of our best selling instruments and kits, should you need some inspiration for your next purchase.

All levels of our product packaging come with the addition of [GS1 Data Matrix labels](#). Switching to the standardised electronic communication system provides complete traceability of our instruments from the point of manufacture to the point of care.

We have worked hard to ensure that all our products and packaging are composed of recyclable materials.



MATERIALS

Instruments

Our instrument handles are Acrylonitrile Butadiene Styrene (ABS) plastic.

The instrument stems and tips are made from Stainless Steel, including our London College Tweezers. Our Examination Tweezers are Steel and Zinc coated.

The coating of our Front Face Mirror is Aluminium, the substrate is Polycarbonate (PC) and the Percussion Tip is Zinc Alloy.

Our Orthodontic Ruler is made from Polypropylene (PP).

In our Complete Periodontal Diagnostic Kits, the Tray is Polypropylene (PP), the Face Masks are a Non-Woven Fabric, the Patient Bib is Paper coated with Polyethylene film, the Tissue Wipe is Paper and the Cotton Rolls are... Cotton.

Sleeves

Our sleeve design comprises of Polyethylene Terephthalate (PET) with a Tyvek® medical dialysing paper strip.

Dispensers and Outer Cartons

The dispenser in which your instruments and kits are stored and the outer carton in which the products are shipped are both made up of 100% recycled paper and cardboard.

Sharps Containers

The Terumo containers we supply are manufactured from high quality, puncture resistant Polypropylene (PP) and Polyethylene (PE).

All the above materials are recyclable.



HOW TO BUY

Online Quoting System

If you would like to obtain an accurate quote inclusive of GST and the delivery charge to your area, our newly launched website features an instant quoting system. All you have to do is visit our [product page](#), add your preferred items to your quote and fill in the required information to receive it instantly.

To make this quicker in future, you can [create an account](#) so you can easily review previous quotes and request new ones without re-entering your details.

Purchasing

When you are ready to place an order, simply contact our team via phone on 07 5445 6486 or [email](#) to purchase items you have been quoted for previously, or to request new ones. Our online quotes are valid for 30 days and can be emailed to our team directly for efficient, accurate order placement.

Samples

If you would like to trial any instruments prior to purchasing, free sterile samples can be requested [here](#).



SHARPS CONTAINERS

As authorised distributors of Terumo sharps containers, we are now offering a clinical waste management service where we can recycle your sharps waste into energy. When you buy a sharps container alongside your instrument purchase, we will deliver the container, collect and recycle your waste for free.

*Waste management service not yet available in all areas of Australia. Please contact our team for more information.

With the aim of simplifying waste collection forecasting for each clinic, the range of containers we have selected to supply hold similar capacities to our instrument dispensers. A single 11L waste container allows for disposal of approximately 300 instruments, the equivalent of 1 instrument dispenser. The 17L waste container allows for disposal of approximately 600 instruments, the equivalent of 2 dispensers.

11L Container with Roll Top Lid
including supply, collection & disposal – \$28.00

17L Container with Roll Top Lid
including supply, collection & disposal – \$35.00

The accessories necessary to install the sharps containers in your clinic are also available. Please [contact](#) our team or your state brand representative for more information.



Single Use Dental Instruments
1/24 Hoopers Road, Kunda Park
QLD Australia 4556
07 5445 6486
sales@sudinstruments.com
www.singleusedentalinstruments.com

