


**SINGLE USE**



**PB92**

**No.12  
& CPITN-C Probe**

## PB92 | No.12 & CPITN-C Probe



### Description

Our sterile combination probes make patient examinations quick and convenient. The ball tip maximises patient comfort when measuring pockets.

### Features

- Combination instrument
- Precision manufactured stainless steel tips
- Ergonomic handle with ridges and grooves for improved grip
- Calibrations; No.12: 3/6/9/12mm; CPITN-C: 3.5/5.5/8.5/11.5mm plus 0.5mm wide ball tip
- Latex-free
- Arrives sterile and ready for use
- GS1 Data matrix
- CE Certified

### Materials

- Stainless steel
- Acrylonitrile butadiene styrene (ABS)